



Product Specification

Balloon Material	Nylon 12
Nominal Pressure	6 atm
Rated Burst Pressure*	16 atm (5 - 6 mm); 15 atm (7 - 8 mm); 14 atm (9 - 10 mm)
Sheath Compatibility	5 Fr (5 - 7 mm); 6 Fr (8 - 10 mm)
Compliance	< 10% between nominal and rated burst pressure
Deflation Time	< 30 seconds
Fold	Tri-fold (\leq 7mm diameter), Penta-fold ($>$ 7mm diameter)
Range of Balloon Sizes	
Diameters	5 to 10 mm
Lengths	20 to 150 mm
Tip Lead in Profile	.048", 1.22 mm
Usable Catheter Length	80 cm \pm 2 cm
Shaft Material	Pebax® 72D, Bi-lumen
Shaft Diameter	0.066", 1.676 mm
Hub	Polycarbonate female y-luer lock
Guidewire	Maximum wire diameter 0.035"
Marker Bands	Platinum/iridium marker bands



*95% Confidence, 99.9% Reliability

Customer is responsible for FDA submission clearance and CE mark prior to product sale.