

**Product Specification**

**Balloon**

Balloon Material	Nylon 12
Nominal Pressure	6 atm
Rated Burst Pressure*	16 atm
Sheath Compatibility	4 Fr
Compliance	< 10% between nominal and rated burst pressure
Deflation Time	< 10 seconds
Fold	Tri-fold

**Range of Balloon Sizes**

Diameters	2 to 4 mm (0.5 mm size increments available)
Lengths	10 to 35 mm

**Catheter Shaft**

Entry Tip Profile	0.019", 0.48 mm
Overall Catheter Length	140 cm ± 2.5 cm usable length
Distal Shaft Material	Nylon 12
Proximal Shaft Material	Pebax® jacketed stainless steel 304
Guidewire	Maximum wire diameter 0.014"


**Hub**

Polycarbonate female luer lock

**Marker Bands**

Swaged platinum/iridium marker bands

\*95% Confidence, 99.9% Reliability



Customer is responsible for FDA submission clearance and CE mark prior to product sale.