


**Balloon Part Number 316048-03020**



Balloon Material	Nylon 12
Nominal Pressure	6 atm
Minimum Burst*	16 atm
Balloon Outer Diameter (A)	3.0 mm
Balloon Length (B)	20 mm
Proximal/Distal Neck I.D. (C/D)	0.036/0.022", 0.91/0.56 mm
Proximal Neck Length (E)	1.5 mm
Distal Neck Length (F)	1.5 mm
Cone Angle Proximal (G)	34°
Cone Angle Distal (H)	34°
Cone Length Proximal (K)	3.4 mm
Cone Length Distal (L)	4.0 mm

\*95% Confidence, 99.9% Reliability  
M316048-03020A