


**Balloon Part Number 316211-05150**



Balloon Material	Nylon 12
Nominal Pressure	6 atm
Minimum Burst*	16 atm
Balloon Outer Diameter (A)	5 mm
Balloon Length (B)	150 mm
Proximal/Distal Neck I.D. (C/D)	0.064/0.050", 1.62/1.27 mm
Proximal Neck Length (E)	5 mm
Distal Neck Length (F)	5 mm
Cone Angle Proximal (G)	40°
Cone Angle Distal (H)	40°
Cone Length Proximal (K)	N/A
Cone Length Distal (L)	N/A

\*95% Confidence, 99.9% Reliability  
M316048-05150A