


Balloon Part Number 316211-06040



Balloon Material	Nylon 12
Nominal Pressure	6 atm
Minimum Burst*	16 atm
Balloon Outer Diameter (A)	6 mm
Balloon Length (B)	40 mm
Proximal/Distal Neck I.D. (C/D)	0.064/0.050", 1.62/1.27 mm
Proximal Neck Length (E)	5 mm
Distal Neck Length (F)	5 mm
Cone Angle Proximal (G)	40°
Cone Angle Distal (H)	40°
Cone Length Proximal (K)	N/A
Cone Length Distal (L)	N/A

*95% Confidence, 99.9% Reliability
316211-06040A

27752 El Lazo Road | Laguna Niguel, CA 92677
 phone 949.448.7056 | fax 949.448.7016 | email info@interfaceusa.com
 ©2009 Interface